# User instruction

The instruction below should be used in conjunction with detailed information on the packaging.

#### Short description of the product

Nitrile protective gloves , powder-free, non-sterile for disposable use

### Full description of the product

run description of the product				
Brand name	: MERCATOR gogrip long			
Reference number	: RP30031002_0001, RP30031003_0001, RP30031004_0001, RP30031005_0001, RP30031006_0001			
Raw material	: nitrile			
External surface	: diamond textured			
Internal surface	: chlorinated			
Cuff	: beaded			
Colour	: Blue			
Shape	: ambidextrous, fitting to the right and left hand			
Size range	: S, M, L, XL, XXL			
AQL	: 1.5			
Quantity in packaging	: 50 pcs. by weight			
Shelf life	: 3 years (from the date of manufacturing)			

#### Storage instructions

is recommended to store the gloves in dry place, in the temperature of 5-35°C and to protect them against direct sunlight and fluorescent light Keep the gloves in a distance of not less than 1m from heating devices, sources of fire and ozone.

Do not keep in direct vicinity of solvents, oils, fuels and lubricants.

#### Food contact

Gloves are marked with food contact symbol  $\overline{X}^0$  and comply with the requirements of European Regulation (EC) No 1935/2004 on materials, European Commission Directive 93/11/EEC and Council of Europe Resolution AP (2004) 4 and with requirement of Regulation (EC) No 2023/2006 on Good Manufacturing Practice.

Result Summary: Test Requested Council of Europe Resolution AP (2004) 4

Extraction conditions (tested for 0.5 hrs in 40°C)	Test Result (limit < 10 mg/dm <sup>2</sup> )
3% Acetic acid	Pass
10% Ethanol	Pass
20% Ethanol	Pass
50% Ethanol	Pass
Rectified Olive oil	Pass

It can be ensured that no substances are transferred to the food which may lead to an alteration of the food. Food disposable gloves are designed for short-term use and frequent changes. It's best to change your gloves regularly.

### PPE classification & compliance

Gloves are category III Personal Protective Equipment as per Annex II of the Regulation 2016/425 and comply to standards: EN ISO 21420: 2020, EN ISO 374-1:2016+ A1:2018 (Type B), EN ISO 374-2:2019, EN16523-1:2015+ A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016

EU Type Examination Certificate issued by: SATRA (Notified Body No. 2777)

Checking of PPE manufactured:



## SATRA Technology Europe Limited

Bracetown Business Park. Clonee.

#### D15YN2P. Republic of Ireland.

Declaration of Conformity and Instruction of Use with the information about the importer are available at: www.mercatormedical.eu

Level1>10min , Level2>30 min, level3>60 min, level4>1 Test results acc. to EN 16523-1:	EN ISO 374-4:2019 Degradation [%]			
Chemical	Level			
40% Sodium Hydroxide (K)	6	-8.1		
37% Formaldehyde (T)	4	14.8		
30% Hydrogen Peroxide (P)	3	31.9		
n-Heptane (J)	4	36.8		

exposure to the challenge chemical. \*\* Permeation rate 7 µg/cm²/min, EN ISO 374-4:2019 degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemica

Test acc. to EN ISO 374-2:2019 – Level 2 (ISO 2859)						
Performance level	Level 3	Level 2	Level 1			
AQL	< 0.65	< 1.5	< 4.0			
Test acc. to EN ISO 374-5:2016						
Prote	ction against bacteria & fungi	Pass				
Protection against viruses		Pass				

EN ISO 374-5:2016 The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimer

#### Intended use

These are non-sterile protective gloves with thicker walls and diamond textured, providing excellent grip and protection when working in oily and wet environment. These gloves are intended for use when there is a need for protection against dangerous substances and infective agents.

#### Precautions and indications for use

Dry hands before putting the gloves on. Before usage, inspect the gloves for any defect or imperfections. These are disposable gloves. Do not let chemical substances get under the gloves through the cuff. If a chemical substance reaches the skin, wash it away immediately with plenty of water with soap. If the gloves get punctured, torn or broken during their use, take them off and put on the new ones. Avoid using gloves diry in the inside as they may cause irritation leading to skin inflammation or more serious damages. The gloves should not be used in contact with open fire and to protect against any sharp tools. The gloves are not intended for welding, electric shock protection, ionizing radiation or from the effect of hot or cold objects.

electric shock protection, jonizing radiation or from the effect of hot or cold objects. This information does not reflect the actual duration of protection in the workplace and the between mixtures and pure chemicals. Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. The chemical penetration resistance has been assessed under laboratory conditions from samples taken from the palm only (except in case where glove is equal to or over 400 mm – where the cuff is tested also) and relates only to the chemical tested and to the tested specimen. It can be different if the chemical is used in a mixture. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test desending on the tamerative abersion and derarchation.

depending on the temperature, abrasion and degradation.

depending on the temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

### Components / hazardous components


Components used in making gloves may cause allergic reactions in some people. Some gloves may contain components known to be a possible cause of allergy for person allergic to them, who may develop contact irritation and/or allergic reaction. In case of an allergic reaction consult a doctor.

#### Disposal


Local regulations regarding the disposal should be applied.

Symbols us	ed on the packaging	5			
PPE	Personal Protective Equipment	NITRILE	Nitrile Gloves		Powder-free gloves
(	For single use only	NON	Non-sterile	$\bigotimes$	Product quality is not ensured if the package is damaged
×	Keep away from sunlight	-	Keep dry	5°C	Temperature limitation 5-35°C
Ŏ <sub>3</sub>	Keep away from ozone	ISO 374-V/TYPE B	EN ISO 374-1 chemical protection gloves; accompanied by a 3-digit code referring to three	150 374-5:2016	EN ISO 374-5 Glove protecting against bacteria, fungi and viruses
ĺ	Consult instructions for use	REF	Catalogue number	LOT	Lot / batch number
#	Model number	${\frown}$	Date of manufacture		Expiry date
	Manufacturer	EC REP	Authorized representative in the European Union		Importer
<b>C€</b> 27	77 CE Mark, PPE	<b>R</b> ï	Suitable for food contact (for details check the Intended use)	PAP	Recyclable packaging
Ì	Package can be treated as municipal waste		Indicates compliance with the requirements of Ukrainian market		Additional information on inner side of package

## How to put the gloves on



# How to take the gloves off



## Manufacturer



## Authorized Representative / Importer

